

PORTABLE ASPIRATORS

Models 405, 300



MODEL 4005



MODEL 300

PORTABLE ASPIRATORS – Gomco® portable aspirators are designed to provide suction that is available when you need it. These units are engineered to supply high performance suction capability at a moment's notice. Model 405 is completely self-contained and has a comfortable handle for easy portability. Model 300 has a Lumitone® enamel finish with an attractive high-impact, fire retardant orange base.

MODEL 405 – This model is designed for general suction use, from emergency room to recovery room, or for outpatient or dental clinic use. It is invaluable for removing mucus from the throats of newborn infants and post-operative patients, for bronchoscopic suction, and for many other bedside uses. The unit is lightweight and compact for convenient storage. The suction system includes

a gauge and suction regulating valve that will create up to 25" Hg vacuum. This unit is UL/ULC listed and includes a hospital grade plug, graduated 1200 ml glass collection bottle or 1500 ml disposable canister, and filter to help prevent fluid and aerosol contamination.

MODEL 300 – This lightweight diaphragm pump unit is designed for general suction use. Its powerful 1/8 horsepower motor creates up to 22" Hg vacuum. Unit includes a gauge and suction regulating valve and is compact for easy storage. This unit is UL listed and includes a hospital grade plug, 600 ml glass bottle or 1100 ml disposable canister, and filter to help prevent fluid and aerosol contamination.

MODEL	COLLECTION BOTTLE	MOTOR	COMPRESSOR	VACUUM RANGE	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
405 4005	Glass Disposable	1/10 hp	Carbon Rotary Vane	0-25 in. Hg 0-635 mm Hg	40 lpm	12 x 9 x 12 in. 30 x 23 x 30 cm	14.5 lbs. 6.6 kg
01-23-0405	Glass	1/20 hp	Dual Diaphragm	0-580 mm Hg	28 lpm	12 x 9 x 12 in. 30 x 23 x 30 cm	14.5 lbs. 6.6 kg
300 3001	Glass Disposable	1/8 hp	Diaphragm	0-22 in. Hg 0-559 mm Hg	29 lpm	12 x 9 x 13 in. 30 x 23 x 33 cm	23 lbs. 10 kg

ORDERING INFORMATION

w/Glass Bottle	w/Disposable Canister	Voltage
01-12-0405	01-12-4005	110-115/50-60
01-23-0405*		230/50-60
01-22-0300	01-22-3001	115/60
01-23-0300		220/50
01-24-0300		220/60

*International only

GOMCO DISPOSABLE COLLECTION UNITS

20-08-0003 Disposable Suction Canister, 1500 ml, w/stem inlet (48/case).
For use with Model 4005 only.